



Slavonska avenija 26/12
Zagreb, Croatia
tel: +385 1 239 68 88
fax: +385 1 239 68 99
e-mail: info@aandb.hr


Biofire® FilmArray® paneli Sindromski pristup u dijagnostici infektivnih bolesti

- **Multiplex PCR paneli**
- Ekstrakcija, amplifikacija, detekcija
- **60** minuta do rezultata
- BioFire testovi za viruse, bakterije, parazite, kvasce i gene antimikrobne rezistencije
- Uzorci: likvor, sinovija, pozitivna hemokultura, BAL, sputum, aspirat, stolica, bris nazofarinksa




 **The BioFire® Respiratory 2.1 (RP2.1) Panel**
Sample Type: Nasopharyngeal swab

 **The BioFire® FilmArray® Gastrointestinal (GI) Panel**
Sample Type: Stool in Cary Blair

 **The BioFire® FilmArray® Pneumonia (PN) Panel**
Sample Type: BAL: (including mini-BAL)
Sputum: (including endotracheal aspirate)

 **The BioFire® Bone and Joint Infection (BJI) Panel**
Sample Type: 0.2 mL of synovial fluid

 **The BioFire® Blood Culture Identification 2 (BCID2) Panel**
Sample Type: Positive Blood Culture

 **The BioFire® FilmArray® Meningitis/Encephalitis (ME) Panel**
Sample Type: Cerebrospinal Fluid (CSF)

Tradicionalno mikrobiološko testiranje



Testiranje sindromskim pristupom PCR metodom

